


SIEMENS

LUMINOS dRF / dRF MAX (90/90-45) SYSTEM TYPICAL ROOM PLAN



The intended use for this Cut Sheet is to communicate the spatial requirements as well as the basic architectural, electrical, structural, and mechanical requirements for this piece of imaging equipment. The information provided in this document is for reference only, during the pre-planning stage, and therefore does not contain any site specific detailed requirements. This information is subject to change without notice. Federal, state and/or local requirements may impact the final placement of the components. It is the customer's responsibility to ensure that the final layout and placement of the equipment complies with all applicable requirements.

LUMINOS dRF / dRF MAX (90/90-45) SYSTEM TYPICAL ROOM PLAN



TYPICAL PLAN

SCALE: 1/8" = 1'-0"

LUMINOS dRF / dRF MAX (90/90-45) SYSTEM SPECIFICATIONS

EQUIPMENT LEGEND								
NO	DESCRIPTION	SMS SYM	WEIGHT (LBS)	BTU/HR TO AIR	DIMENSIONS (INCHES)			REMARKS
					W	D	H	
①	CONTROL CONSOLE (GENERATOR & REMOTE CONTROL)	Ⓢ	10	34	19 1/2	13 3/4	3 7/8	ON CUSTOMER'S COUNTER
②	FLUOROSPOT COMPACT – KEYBOARD AND MOUSE	Ⓢ	–	–	–	–	–	ON CUSTOMER'S COUNTER
③	B/W FLAT SCREEN CONTROL ROOM MONITOR	Ⓢ	–	–	–	–	–	ON CUSTOMER'S COUNTER
④	FLUOROSPOT-COMPACT CONTAINER (UNDER COUNTER)	Ⓢ	110	1,468	17 3/8*	32 1/2*	27*	*INCLUDES 5 1/8" BEHIND, 4" ABOVE AND ALONG CONTAINER.
⑤	POLYDOROS F80 (80 kW) GENERATOR CABINET	Ⓢ	838	2,048**	31 1/2	17 1/8	86 3/4	**DURING OPERATION, 1,195 IN STANDBY MODE
⑥	LUMINOS dRF (+90/-90) REMOTE TABLE	Ⓢ	2,911	2,730	83	75	107***	***MAX. AT 60" SID.
⑦	18" FLAT DISPLAY ON CART	Ⓢ	187	256	33 1/4	28 3/8	60	
⑧	YSIO WALL STAND WITH MOBILE DETECTOR (LEFT LOADING)	Ⓢ	551	819	30	37*A	83	*A – MAX. IN HORIZONTAL POSITION
⑨	(*) DOCKING STATION (WALL MOUNTED)	Ⓢ	40	256	20 3/8	7 13/16	16 1/4	WITHIN 11.5 FT. OF WALL STAND
⑩	4.25M RAILS FOR FULLY SYNCHRONIZED TUBE STAND	Ⓢ	59	–	167 3/8	3	3 1/2	SIZE AND WEIGHT PER RAIL
⑪	3M TRANSVERSE BRIDGE AND X-RAY TUBE STAND	Ⓢ	772	853	119 1/4	39	*43	*TUBE CARRIAGE IN PK POS
⑫	GRID HOLDER (WALL MOUNTED)	Ⓢ	22	–	21 11/16	4	16 9/16	SUGGESTED LOCATION
⑬	(MAX) DVI INTERFACE FOR EXAM ROOM DISPLAY (OPTION)	Ⓢ	–	–	–	–	–	CUSTOMER PROVIDES MONITOR
⑭	(MAX) CHARGING STATION FOR MAX DETECTORS (OPTION)	Ⓢ	–	–	–	–	–	ON CUSTOMER'S COUNTER
⑮	(MAX) ACCESS POINT (TOP OF GENERATOR)	Ⓢ	4	–	8	5.5*A	3*B	*A 12 1/4, *B 9" INCLUDING ANTENNAE
Ⓢ	(*) – THIS ITEM IS ELIMINATED FOR ALL MAX SYSTEMS (MAX) – ITEM ADDED FOR MAX SYSTEMS ONLY		–	–	–	–	–	

CEILING HEIGHT REQUIREMENTS		
SYSTEM CONFIGURATION		CEIL. HT. RANGE
TABLE ONLY		8'-3" TO 10'-6"
CEILING STAND + TABLE +/- WALL STAND		8'-9" TO 9'-3"
CEILING STAND + TABLE +/- WALL STAND (+TUBE EXT.)		9'-4" TO 10'-2"
PATIENT TABLE (ONLY) EXAMS		MINIMUM CEIL. HT.
TABLE X-RAY TUBE TO INTEGRATED TABLE BUCKY (0 DEG. POS.) USING LOWEST POSSIBLE TABLE TOP HEIGHT 2'-1 1/2" [25 1/2"]	45" SID	8'-3"
	60" SID	
TABLE X-RAY TUBE TO UPRIGHT (+/-90) INTEGRATED TABLE BUCKY FOR ORTHO (NOTE 1)		8'-10"
1. WITH 8'-3" MIN. CEILING HT., +/-90 DEG. VERTICAL TABLE TOP TRAVEL IS 13".		

LUMINOS dRF TECHNICAL DATA	
ENVIRONMENT:	59°F – 95°F OPERATING ROOM TEMPERATURE 20% – 75% PERMISSIBLE RELATIVE AIR HUMIDITY (NON-CONDENSING)
TRANSPORTING/RIGGING:	LARGEST CRATE: 97"L X 60"W X 56"H MINIMUM ELEVATOR SIZE: 117"L X 34"W X 60"H HEAVIEST SINGLE PIECE: 2,448 LBS. WITH PACKING 1,874 LBS. WITHOUT PACKING MINIMUM DOOR OPENING (FOR TABLE): 48" WIDE WITH MINIMUM 6'-11" CORRIDOR WIDTH. ALSO, 34" WIDE OPENING WITH MINIMUM 9'-7" CORRIDOR WIDTH. LARGEST PIECE (TABLE TRANSPORT CARRIAGE): 115"L X 33"W X 57"H (WHEELS OUTSIDE) 106"L X 33"W X 57"H (2-WHEELS INSIDE) 94"L X 33"W X 57"H (ALL WHEELS INSIDE)
REV 0	

MINIMUM CEILING HEIGHT	CEILING HEIGHT WITHOUT RESTRICTION	RECOMMENDED CEILING HEIGHT
W/RESTRICTION	RESTRICTION	
SEE CHART	10'-7"	9'-5"

LUMINOS dRF / dRF MAX (90/90-45) SYSTEM SPECIFICATIONS

POLYDOROS F80 80kW	
X-RAY GENERATOR POWER REQUIREMENTS	
INCOMING POWER:	480 VOLTS, 3 PHASE, 60Hz
CIRCUIT BREAKER:	80 AMPS.
GENERATOR OUTPUT:	80 kW
ALLOWABLE IMPEDANCE:	≤ 0.16 Ω
MAXIMUM MOMENTARY LOAD:	126 kVA
LINE VOLTAGE VARIATION:	± 10% MAX.
PHASE IMBALANCE:	± 2%
FREQUENCY VARIATION:	± 1 Hz
<p>NOTE: ALL INCOMING POWER SUPPLIES, FOR THE SIEMENS EQUIPMENT, ARE TO BE DEDICATED (BACK TO SOURCE) ISOLATED AND INSULATED FROM ANY OTHER EQUIPMENT, SUCH AS, ELEVATORS, GENERATORS, HVAC SYSTEMS, ETC.</p> <p>A NEUTRAL CONDUCTOR, IF PRESENT, IS NOT USED FOR THE LINE VOLTAGE CONNECTION TO THE SIEMENS EQUIPMENT. IF THE NEUTRAL CONDUCTOR IS PROVIDED, IT SHOULD NOT BE ELECTRICALLY CONNECTED AT ANY POINT IN THE POWER DISTRIBUTION TO THE SIEMENS EQUIPMENT UNLESS SPECIFICALLY REQUIRED. UNINTENTIONAL NEUTRAL TO GROUND BONDS MAY VIOLATE LOCAL AND NATIONAL ELECTRICAL CODES, AS WELL AS CREATE GROUNDING PROBLEMS.</p> <p>IF AN ON-SITE TRANSFORMER IS REQUIRED TO OBTAIN XP MODALITY OPERATING VOLTAGE, IT MUST BE OF SUFFICIENT CAPACITY AND CHARACTERISTICS TO MAINTAIN SUPPLY VOLTAGE AND IMPEDANCE REQUIREMENTS (TRANSFORMER & CONDUCTORS).</p>	
<p>ATTENTION: SIEMENS MEDICAL SYSTEMS, INC. RECOMMENDS THAT THE INCOMING POWER LINES BE ANALYZED WITH RESPECT TO TRANSIENT SURGES AND IMPULSES, SAGS, AND OVERVOLTAGES.</p>	
REV 2	

WIRELESS DETECTOR CONNECTION
<p>OPERATION OF THE WIRELESS DETECTOR CAN BE AFFECTED BY OTHER WLAN DEVICES IN THE VICINITY OF THIS INSTALLATION. TO AVOID ANY CONFLICTS, THE CUSTOMER MUST PROVIDE A LIST OF EXISTING WLAN CHANNELS (FREQUENCIES) OR THE SPECIFIC CHANNEL (FREQUENCY) THEY DESIRE TO BE USED FOR THE WIRELESS DETECTOR.</p>
<p>THE WIRELESS CONNECTION IS ENCRYPTED (WPA2) AND IS BASED ON TWO WLAN STANDARDS, WITHIN WHICH SEVERAL CHANNELS (FREQUENCIES) ARE AVAILABLE:</p> <ol style="list-style-type: none"> 1) 11G STANDARD – OPERATES AT 2.5 GHz 2) 11A STANDARD – OPERATES AT 5 AND 6 GHz
<p>THE STANDARD (11G OR 11A) CAN BE SET BY SIEMENS SERVICE VIA THE SERVICE SOFTWARE INSTALLED ON THE IMAGING SYSTEM.</p>
<p>THE WIRELESS CONNECTION IS ONLY USED TO TRANSFER DATA BETWEEN SIEMENS EQUIPMENT AND IS NOT USED TO SEND DATA TO THE CUSTOMER'S NETWORK.</p>
REV 0

SYSTEM TECHNICAL DATA			
TRANSPORTING INFORMATION			
		SIZE	WEIGHT
TRANSVERSE BRIDGE	3 M	126"L x 32"W x 10"H	419#
	4 M	174"L x 32"W x 10"H	512#
LONGITUDINAL RAILS	4 M	167"L x 3"W x 4"H	59# EACH
	5 M	197"L x 3"W x 4"H	82# EACH
DCS-1/2	4 M	167"L x 3"W x 3"H	441#
X-RAY TUBE STAND (FULLY SYNCHRONIZED)		67"L x 41"W x 53"H	827#
WALL STAND WITH MOBILE DETECTOR (WITH PACKING)		35"L x 93"W x 41"H	898#
MINIMUM DOOR OPENING:	3'-5 3/8"		
MINIMUM CORRIDOR WIDTH:	6'-11"		
ENVIRONMENTAL CONDITIONS			
	IN OPERATION	TRANSPORT	
PERMISSIBLE AMBIENT TEMPERATURE (WITH WIRELESS DETECTOR)	59°F TO 82°F	-4°F TO 131°F	
PERMISSIBLE RELATIVE HUMIDITY	20% TO 75%	5% TO 95%	
		REV 0	

REMOTE SYSTEM DIAGNOSTICS
<p>SIEMENS REMOTE SERVICES (SRS) REQUIRES A CONNECTION BETWEEN THE SRS REMOTE SERVER AND SIEMENS SYSTEMS VIA REMOTE LOCAL AREA NETWORK ACCESS, TO ENSURE THE UPTIME OF YOUR SYSTEM. A CUSTOMER VPN CAPABLE FIREWALL OR OTHER VPN APPLIANCE IS PREFERRED.</p>

FOR MORE INFORMATION
<p>FOR MORE DETAILED PLANNING REQUIREMENTS FOR THIS SYSTEM, SEE THE TYPICAL FINAL DRAWING SET NUMBER: 08015</p>

LUMINOS dRF / dRF MAX (90/90-45) SYSTEM SPECIFICATIONS

MAXIMUM CABLE DISTANCES BETWEEN COMPONENTS								
	CONTROL CONSOLES	FLUOROSPOT COMPACT	GENERATOR	LUMINOS dRF TABLE	MOBILE FLAT DISPLAY CART	CEILING TUBE STAND	DETECTOR WALL STAND	DOCKING STATION
GENERATOR	-	59'-0"	-	22'-0"	-	32'-0"	36'-0" 52'-6" (NOTE 1)	36'-0"
FLUOROSPOT COMPACT	11'-0"	-	59'-0"	-	32'-0"	-	-	-
LUMINOS dRF TABLE	59'-0"	-	22'-0"	-	-	-	-	-

NOTE 1 - MAX SYSTEM OPTION: 16 METER CABLE EXTENSION MUST BE SELECTED.

THE DISTANCES LISTED ABOVE ARE CALCULATED AS THE MAXIMUM CABLE LENGTH BETWEEN CABLE ENTRY POINTS. DEPENDING ON THE COMPONENT, THE CABLE ENTRY POINT MAY BE IN FLOOR, WALL OR CEILING. VARIOUS ARRANGEMENTS OF COMPONENTS ARE POSSIBLE AS LONG AS THE DISTANCES SHOWN ARE MAINTAINED AND THE SYSTEM FUNCTIONALITY IS NOT ADVERSELY AFFECTED.